

CHOCOLATE FOR TWO BASKET

1. Crumple two pieces of tissue on to the tray.
2. From left to right - place the acetate of assorted chocolates, then the seafoam and chocolate almonds, skip a space as wide as the ribbon or fill it in with another piece of tissue latter – last add the acetate sampler to the right side.
3. Make a 30inch long bag and gather it at the top securing it with wire. Shrink the bag around the product only being careful not to heat the top knot.
4. Add a big bow and loose tails.
5. The seafoam and the chocolate covered almonds can be changed according to customer preference but the two variety packs should always remain.